

## Cardboard Freezer Boxes

### I. APPLICATIONS

1. Storage and organization of any range of 0.5ml to 50ml tubes
2. Great for upright or chest freezers

### II. FEATURES

1. Made from durable cardboard with variety of outer coatings: (Superior white , Plasti-coat, Cryo-color)
2. Chemically resistant to alcohols and mild organic solvents
3. Temperature range: stable from -196°C to 121°C
4. Suitable for liquid nitrogen storage

### III. BOX DESIGN

1. Available in 2 inch and 3 inch sizes with either 81-well, 100-well, 36-well or 16-well dividers
2. Select from our Superior white coating, our Plasti-coat laminated coating, our Cryo-color or ID-color assorted colored boxes. (Pink, blue, yellow and green)
3. Holds 0.5ml, 1.5ml or 2.0ml microtubes, and 5ml, 15ml or 50ml large volume tubes
4. Removable lid includes a inventory list writing area
5. Cardboard dividers

### IV. PACKAGING

1. Bulk, 100 boxes/case

### V. DESCRIPTION

Item	Description	Packaging	Sterile
90-1536	36-well, superior white	100 pieces/case	No
90-5016	16-well, superior white	100 pieces/case	No
90-1281	2in, 81-well, no outer coating	100 pieces/case	No
90-1200	2in, 100-well, no outer coating	100 pieces/case	No
90-2281	2in, 81-well, plasti-coat	100 pieces/case	No
90-2200	2in, 100-well, plasti-coat	100 pieces/case	No
90-2381	3in, 81-well, plasti-coat	100 pieces/case	No
90-2300	3in, 100-well, plasti-coat	100 pieces/case	No
90-5281	2in, 81-well cryo-color, assorted	100 pieces/case	No
90-5200	2in, 100-well cryo-color, assorted	100 pieces/case	No

### VI. STERILITY

1. N/A


## VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

## VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


## IX. PRODUCT PHOTOS



**X. PRODUCT PARAMETER**





Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2381	Box	133.0	125.0	15.0	15.0	71.0	69.2	55.0
		<b>H</b>	<b>I</b>	<b>J</b>	<b>K</b>	<b>L</b>	-	-
		52.5	51.0	133.0	130.0	76.0	-	-







Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2300	Box	133.0	125.0	12.5	12.5	71.0	69.2	55.0
		<b>H</b>	<b>I</b>	<b>J</b>	<b>K</b>	<b>L</b>	-	-
		52.5	51.0	133.0	130.0	76.0	-	-







**Dimensions=mm**

Cat #	Description	A	B	C	D	E	F	G
90-1281	Box	133.0	127.0	15.0	15.0	45.5	40.0	30.0
90-2281		H	I	J	K	L	-	-
90-5281		34.0	32.5	133.0	130.0	53.0	-	-







**Dimensions=mm**

Cat #	Description	A	B	C	D	E	F	G
90-1200	Box	133.0	127.0	12.5	12.5	45.5	40.0	30.0
90-2200		H	I	J	K	L	-	-
90-5200		34.0	32.5	133.0	130.0	53.0	-	-



**Dimensions=mm**


Cat #	Description	A	B	C	D	E	F	G
90-1536	Box	145.0	142.0	25.0	25.0	75.0	73.5	70.0
		<b>H</b>	<b>I</b>	<b>J</b>	<b>K</b>	<b>L</b>	-	-
		68.5	67.5	148.0	145.0	76.0	-	-



90-1536

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-5016	Box	145.0	142.0	35.0	35.0	75.0	73.5	70.0
		<b>H</b>	<b>I</b>	<b>J</b>	<b>K</b>	<b>L</b>	-	-
		68.5	67.5	148.0	145.0	76.0	-	-



90-5016