

Product Name:	Centrifuge Tube
Product Application:	Googlab centrifuge tubes could be commonly used to collect cells and separate biomacromolecules by centrifugation in biochemistry and other biological fields. They can also be used to store photosensitive samples and special chemiluminescence reagent.
Product Materials:	PP (tube) + PP(cap) Meets USP, Class VI standards
Product Specifications:	Temperature range: stored at room temperature Shelf life: 3 years after date of production. (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137.
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Sealing test: qualified Packaging inspection: qualified
Features:	<ol style="list-style-type: none"> 1. Printed clear graduations with marking area. 2. Polypropylene tube and high density polyethylene cap. 3. High centrifugal strength up to 12,000xg. 4. Withstand temperature range: Approx.-80°C to 121°C. 5. Flat markable cap with leakproof seal. 6. One-hand operation design for the caps is easy to handle. 7. Sterilized by E-beam, SAL=10⁻⁶. 8. Non-Pyrogenic, DNase/Rnase free. 9. Tubes are available in two usage levels: Premium Level and Standard Level.

Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

Product Range:

Centrifuge Tube


Cat. No.	Volume (ml)	Color	Max.RCF (xg)	Temperature Range	Type of Cap	Pack	Case
G66010522	15	Transparency	12000xg	-80°C do +121°C	Flat Top	10x50pcs.	500pcs.
G66020522	50	Transparency	12000xg	-80°C do +121°C	Flat Top	10x50pcs.	500pcs.

Technical Drawing:

Centrifuge tube

G66010522


• 15ml



Centrifuge tube

G66020522

• 50ml



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com