

<b>Product Name:</b>	<b>PCR Rack</b>
<b>Product Application:</b>	It is used to place the PCR tubes, which is convenient for adding liquid to the tubes and keeps the PCR tubes in a vertical position. Fits all NEST PCR tubes.
<b>Product Materials:</b>	polypropylene (PP)
<b>product Specifications:</b>	Temperature range: stored at room temperature Shelf life: 5 Years after date of production ( ensure package is in good )
<b>Sterilization:</b>	No
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: pass <span style="float: right;">Packaging inspection: pass</span>
<b>Features</b>	<ol style="list-style-type: none"> <li>1. Made from polypropylene</li> <li>2. Fit for fixing all NEST PCR tubes.</li> <li>3. Stackable</li> <li>4. Transparent polycarbonate lid for well identification</li> <li>5. Blur Rack &amp; Transparent Lid</li> <li>6. Withstand temperature range: Approx.-80°C to 121°C</li> <li>7. Autoclavable</li> </ol>

**Product Range**

Cat.No.	Spec	Cap Type	/Pack	/Case
407001	8*12	Flat	1 Rack & 1 Lid	25 Rack & 25 Lid
407101	8*12	Flat	5 Racks & 1 Lid	25 Rack & 5 Lid



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

## Product picture



## Technical Drawing



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)