

Product Name:	Pasteur Pipette	
Product Application:	It is suitable for drawing, transferring and blending a certain volume of liquid. Both sterile and non-sterile Pasteur pipettes are available.	
Product Materials:	PP	
product Specifications:	Temperature range: stored at room temperature Shelf life: 2 Years after date of production (ensure package is in good)	
Sterilization:	Yes.	
Pyrogens:	Non-Pyrogenic	
RNase/DNase:	DNase/RNase free	
BSE/TSE:	These products are deemed animal free.	
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: pass Packaging inspection: pass	
Features	<ol style="list-style-type: none"> NEST Pasteur pipettes are made of PP and manufactured in ISO 8 grade cleaning room. Good chemical stability. Transparent with clear graduations for easy observation The tube wall has good fluidity, good elasticity and easy control. It is suitable for rapid transfer of liquids. Pipette can be twisted for use with micro vessels and special-shaped containers. High precision and good repeatability. Individually packed pasteur pipettes are sterilized by EO. DNase/RNase, endotoxin and pyrogenic free. 	

Product Range


Cat. No.	Volume(mL)	Length(mm)	Pack	Sterile	/Pack	/Box	/Case
318012	1	147	individually packed	YES	1	500	2000
318112	2	155		YES	1	500	2000
318212	3	156		YES	1	500	2000
318314	3(Extra Long)	177		YES	1	200	2000
318513	5	200		YES	1	150	1500
318415	10	285		YES	1	100	1000
318031	1	147	Bagged	NO	/	500	5000
318131	2	155		NO	/	500	5000
318231	3	156		NO	/	500	5000

Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com


Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Technical Drawing


1ml




2ml




3ml




3ml (Extra Long)



5ml




10ml



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Product Picture



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)