



Product Name:	Serological Pipette
Product Application:	Googlab serological pipettes can be used for liquid transfer in biological experiments such as cell culture and bacterial culture.
Product Materials:	Polystyrene (PS).
Product Specifications:	Temperature range: stored at room temperature Shelf life: 3 years after date of production
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137.
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Stress test: qualified Sealing test: qualified Packaging inspection: qualified
Features:	<ol style="list-style-type: none"> 1. Made from high clarity polystyrene. 2. Withstand temperature range: Approx. 0 to 40 2. Individual packaging. 3. Edged color-coded band for identifying Pipettes of different volumes. 4. Negative graduations for extra volume. 5. Sterilized by E-beam. 6. Non-pyrogenic, non-cytotoxic and non-hemolytic. 7. DNase/RNase free. 8. Polyolefin fiber filter to reduce contamination.

Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com


Product Range:

Serological pipettes


SKU	Volume (ml)	Sterile	Scale	LENGTH	Type of Packaging	Pack	Case
G33240011	1	Yes	Yes	280mm	Individual	6x100pcs.	600pcs.
G33250011	2	Yes	Yes	280mm	Individual	6x100pcs.	600pcs.
G33260011	5	Yes	Yes	343mm	Individual	4x50pcs.	200pcs.
G33270011	10	Yes	Yes	343mm	Individual	4x50pcs.	200pcs.
G33280011	25	Yes	Yes	312mm	Individual	4x50pcs.	200pcs.
G33290011	50	Yes	Yes	350mm	Individual	5x25pcs.	125pcs.

Technical Drawing:


1ml Serological pipettes



2ml Serological pipettes



5ml Serological pipettes



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland


E-mail: office@genoplast.com

Warehouse in the USA


GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com


10ml Serological pipettes



25ml Serological pipettes



50ml Serological pipettes



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com